

Any Label. Any Printer. Your Design

1-Click Print[®]

Faster, Smarter Printing



NEW Auto Fill Features:

- Manufacturer Lot No.
- Manufacturer Expiration Date
- GTIN
- Bottle Serial no.
- Beyond Use Date



No Manual Data Entry Needed

Package/Label Step 1 - Please enter the Product ID.
Enter Product ID [0004006832] [Clear]

Step 2 - Please enter the packaging information.


Auto Enters: Mfg Lot # (Manufacturer Lot #)
Mfg Exp. Date (Manufacturer Exp Date)
Auto Calculates BUD! (Beyond Use Date)

1-Click Print[®]

Auto Enters: Bottle Qty (Package/Label Quantity)

Previous Label Recall
Hazardous Drug Warning
Drug Image
Control Drug Symbol

Print Package/Label [Print] [Print Preview] [Close]



The software that powers iPack[®] Rx

clearTag[®]

Bar Code Labels

- Clear stemmed flag labels for vials, ampules, and syringes
- Scratch and moisture resistant
- Clear stem will not obscure manufacturer's labeling



m:Print[®] Features:

Custom Label Templates:

Design unlimited custom templates for roll, sheet, and flag labels with a few clicks of a mouse. Works with all label supplies.

Any Label. Any Printer. Any Size. Your Design.

Design custom labels and m:Print[®] will print it with any Windows[®] supported printer.

Automatic Job Numbering:

Software assigns unique identifier for each packaging and labeling run. No manual logs needed.

Third Party Drug Database:

NDC drug database with over 150,000 drug records, including complete drug information, drug images, and tallMAN lettering for look and sound-alike drugs.

Bar Code Constructor:

Automatically build DSCA Bar Codes with a check box. Build Custom 2-D and Linear Bar Codes with multiple drug data fields within the bar code.

2-D and Linear Bar Codes:

Data Matrix, Aztec, and most linear bar codes.

Network Version Option

Multiple Report Options:

Advanced data filtering of dates and fields.

Microsoft SQL Server:

Robust enterprise-class database makes sure your data stays secure and safe.

Minimum Computer Specifications:

- Windows[®] 10 Pro, 64 Bit
- 8 GB RAM, 128 GB Hard Drive

1-Click Print[®]

Faster, Smarter Printing

with

Rules-based Workflow

making error-free medication labeling easier than ever

New in v4

- Scan and Print DSCSA Serialized Bar Codes
- Automatic Entry of Mfgr Exp Date, Lot #, and Serial #
- Automatic Calculation of BUD
- USP 1146 & USP 1178 BUD Calculations
- Automatic Previous Label Template Recall
- Automatic Package Numbering (w/ Repeat # option)
- Automatic Job Numbering
- Improved Bar Code Constructor
- Automatically Build DSCA Bar Codes
- Build Custom 2-D and Linear Bar Codes
- Hazardous Drug Warning
- Network Version Option
- Remote RPh Approval
- Enhanced Reports
- Multiple Preformatted Reports
- Pill & Label Image printed on reports
- Pre-Print RPh Approval
- Enhanced Data Filtering Functions
- Welcome Screen with User Message Board
- More Customizable Data Fields
- Scheduled Database Backup
- Import Custom Drug Information
- Fractional Pill Indicator
- Multi-Language Enabled for International Use



Pearson Medical
Technologies SAFE • SMART • EFFICIENT



2804 North Bolton Avenue
Alexandria, LA 71303
Phone: (318) 619- 1129
Fax: (318) 473- 2879
info@pearsonmedical.com

Phone 866-640-3603
www.pearsonmedical.com